

# Fantastic Voyage: A new nanoscale view of the biological world

**Echoing the journey through the human body in *Fantastic Voyage*, doctors might soon be able to track individual donor cells after a transplant, or to find where and how much of a cancer treatment drug there is within a cell.**

New technology described in a study published today in the open access journal *Journal of Biology* makes it possible to image and quantify molecules within individual mammalian or bacterial cells. Claude Lechene and colleagues describe the development of multi-isotope imaging mass spectrometry (MIMS), which has applications in all fields of biology and biomedical research.

"This method allows us to see what has never been seen before, and to measure what has never before been measured," Lechene says. "Imagine looking into a building, slice by slice. You can see not only that it contains apartments, but also that each apartment contains a refrigerator. You can see that there are tomatoes in the refrigerator of one apartment, and potatoes in the refrigerator of another. You can count how many there are and measure how fast they are used and replaced. It is this level of resolution and quantification that MIMS makes possible within cells."

Lechene, of Harvard Medical School and Brigham and Women's Hospital in the US, worked with colleagues from around the world to develop and test the new methodology.

A beam of ions is used to bombard the surface atoms of the biological sample, and a fraction of the atoms are emitted and ionized. These "secondary ions" can then be manipulated with ion optics – in the way lenses and prisms manipulate visible light - to create an atomic mass image of the sample. Lechene et al. developed MIMS by combining the use of a novel secondary-ion mass spectrometer developed by Georges Slodzian, from the Universit#1081; Paris-Sud in France, labeling with stable isotopes and building quantitative image-analysis software.

MIMS can generate quantitative, three-dimensional images of proteins, DNA, RNA, sugar and fatty acids at a subcellular level in tissue sections or cells. "Using MIMS, we can image and quantify the fate of these molecules when they go into cells, where they go, and how quickly they are replaced," says Lechene.

The method does not need staining or use of radioactive labelling. Instead, it is possible to use stable isotopes to track molecules. For example, researchers could track stem cells by labelling DNA with <sup>15</sup>N. "These stable isotopes do not alter the DNA and are not toxic to people; with MIMS and stable isotope labelling we could track these cells, where they are and how they have changed several years later," says Lechene.

"The most significant feature of this technique is that it opens up a whole new world of imaging; we haven't yet imagined all that we can do with it," says Peter Gillespie from the Oregon Health and Science University in Portland, USA in an accompanying news article, also published today in *Journal of Biology*.

*Citation:* High-resolution quantitative imaging of mammalian and bacterial cells using stable isotope mass spectrometry Claude Lechene, Francois Hillion, Greg McMahon, Douglas Benson, Alan M Kleinfeld, J Patrick Kampf, Daniel Distel, Yvette Luyten, Joseph Bonventre, Dirk Hentschel, Kwon Moo Park, Susumu Ito, Martin Schwartz, Gilles Benichou and Georges Slodzian *Journal of Biology* 2006, 5:20 (5 October 2006)

Source: BioMed Central

*This document is subject to copyright. Apart from any fair dealing for the purpose of private study, research, no part may be reproduced without the written permission. The content is provided for information purposes only.*